

Department of Chemistry

Name: Kavlika Singh

Class Roll: 21

[illegible]

① Name:- Kartika Singh Roll no:-21 Date:- 5-10-23

Determination of chloride content of water:-



Aim:-
Determination of chloride ion in a given water sample by Argentometric method (using chromate indicator solution)

Theory:-

When Silver nitrate solution is added to a given solution of chloride containing few drops of indicator (K_2CrO_4), white silver chloride is precipitated initially

[Since $K_{sp}(AgCl) < K_{sp}(Ag_2CrO_4)$, K_{sp} : Solubility product].

After finishing all the chloride ion as $AgCl$, Ag^+ binds CrO_4^{2-} ion and forms a brick red precipitate of Ag_2CrO_4 . This brick red colour indicates the end point.



Apparatus:-

- Burette
- Pipette (25ml)
- Conical flask

Reagents:-

- $AgNO_3$ solution (N/50)
- K_2CrO_4 indicator

①

Determination of Water Quality Parameter (Red-ox Titration)

AIM

Estimation of Iron present in Mohr's salt solution using standardized KMnO_4 solution.

THEORY

In a red-ox titration, a substance is oxidized and the other is reduced. In other words, Oxidation & reduction occurs simultaneously. Oxidation is the process of loss of one or more electrons and reduction is the gain of electrons by atoms or ions. The reagent undergoing reduction is called oxidizing agent (oxidant); and the reagent undergoing oxidation is called reducing agent (reductant). In redox titration, a reducing agent is titrated against an oxidizing agent and vice-versa. Normally potassium permanganate and iodine are commonly used as oxidizing agents and maleic acid, Ferrous ammonium sulphate (Mohr's salt), sodium thiosulphate, sodium oxalate etc. are commonly used reducing agents. Potassium permanganate is a powerful and commonly available oxidizing agent.

④

Discussion:-


1. After addition of silver nitrate from burette into the solution, silver chloride is precipitated.
2. Solubility product of silver chromate is greater than silver chloride, thus silver chloride is precipitated until all the chloride ions are attached to the silver forming silver chloride.
3. The brick red colouration of the solution is due to formation of silver chromate.

Conclusion:-

The amount of chloride ion (determined by Argentometric method) presents in given water sample is 17.04 mg/lit.

Date of Experiment:- 05-10-23.

Date of Submission:- 27-11-23

Signature of Teacher-  27/11/23

Results and Calculations:-

Table : Determination of Chloride Ion:-

SL No.	Volume of Sample taken (ml)	Burette reading		Volume of AgNO_3 (ml)	Concordant volume of AgNO_3 (ml)	Amount of Chloride ion present (PPM)
		Initial	Final			
1	25	45.0	45.6	0.6		
2	25	45.6	46.2	0.6	0.6	17.04
3	25	46.2	46.9	0.7		mg/lit

Calculation

$$\text{A.C.B. } V_1 S_1 = V_2 S_2$$

where, V_1 = volume of water ~~sample~~
Sample taken (25 ml)

$$\rightarrow 25 \times S_1 = 0.6 \times \frac{1}{50}$$

$$25 \times S_1 = 0.012$$

$$S_1 = 0.00048$$

Hence, the strength of water sample due to presence of Cl^- is 0.00048.

S_1 = Strength of water sample (due to presence of Cl^- ion = ?)

V_2 = volume of AgNO_3 solⁿ (0.6)

S_2 = Strength of AgNO_3 solⁿ
 $\text{N}/50$ or $1/50$.

1000 ml 1(N) water sample due to presence of Cl^- ions contains 35.5 gm of Cl^- ion.

$$1000 \text{ ml of } 1 \text{ N water sample} = \frac{0.00048 \text{ (N)}}{0.00048 \times 35.5 \times 1000} = 17.04 \text{ mg/lit}$$

Procedure :-

1. The burette was filled with AgNO_3 (N/50) solution.
2. The pipette was washed with distilled water first and then rinsed with given sample of tap water. 25 ml of sample of tap water was pipetted out in a clean conical flask.
3. 2 drops of K_2CrO_4 indicator were added.
4. AgNO_3 was added from the burette, shaking the flask constantly. A white precipitate of AgCl was obtained. After the addition of more drops of AgNO_3 solution, a red colour appeared in the flask as local point which disappeared upon shaking.
5. The addition of AgNO_3 solution was continued till solution turned reddish brown (brick red colour). The readings were taken.

we know,

$$V_1 S_1 = V_2 S_2$$

$$10.1 \times 0.098 = 10 \times S_2$$

$$S_2 = 0.09898 \text{ N}$$

1000 cc 1(N) Mohr's salt solution contains 56 gm of Fe.

1000 c.c 0.09898(N) Mohr's salt solution contains
= 56×0.09898 gm of Fe

$$= 5.5428 \text{ g (L)}$$

Discussion:-

1. Oxidation reduction takes place simultaneously in this experiment.
2. KMnO_4 solution is standardized with the help of standard Oxalic acid using the formula $V_1 S_1 = V_2 S_2$
3. KMnO_4 acts as a self-indicator in this experiment.
4. KMnO_4 acts as a strong oxidizing agent in presence of H_2SO_4 thus H_2SO_4 is added in the solution.
5. Strength of Mohr's salt solution is determined with standardized KMnO_4 solution using formula $V_1 S_1 = V_2 S_2$
6. When all the reducing agent has been oxidized then the excess drop of KMnO_4 resulted in change of the colour.

Result & Calculation :-

Table :- Titration of Mohr's salt solution Against Standardized KMnO_4 solution (Redox - Titration) :-

No of observation	Vol. of Mohr's salt solution (ml)	Burette Reading		Vol. of KMnO_4 sol ⁿ consumed (ml)	concordant reading (ml)	Strength of KMnO_4 solution (N)	Strength of Mohr's salt solution (N)	Iron present (g/100)
		Initial	Final					
1	10	0.0	10.1	10.1	10.1	0.098	0.098	5.54%
2		10.1	20.2	10.1				
3		20.2	30.4	10.2				

y. Datta
12.10.23

From the normality equation, $V_1 S_1 = V_2 S_2$ the value of S_2 is determined.

Where, V_1 = volume of KMnO_4 solution = 10.1 ml

S_1 = Strength of KMnO_4 solution = 0.098 N

V_2 = volume of Mohr's salt solution = 10 ml

S_2 = Strength of Mohr's salt solution

(3)

with the strength of the Mohr's salt solution.

Apparatus:-




- Graduated pipette (10 ml)
- Burette
- Conical Flask (250 ml)
- Measuring cylinder (10 ml)

Reagents:-


- Mohr's salt solution [$\text{FeSO}_4, (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$]
- Standardized Potassium permanganate solution
- 2(N) Sulfuric acid

(2)

and its oxidizing power depends on the acidity of the solution. In most of these titrations, solution of reducing agents should be made acidic before carrying out the titration. In acidic solution, reduction of potassium permanganate is represented by :-



so the complete equation is



Procedure :-

Known strength of KMnO_4 solution, 2(N) H_2SO_4 solution and unknown strength of Mohr's salt solution are supplied. Now the burette is rinsed with KMnO_4 solution & then filled with it. After that, 10 ml of Mohr's salt solution was pipetted out in a 250 ml conical flask & 10 ml of dilute H_2SO_4 solution was added to it and then it was titrated by standardized KMnO_4 solution. At the end point, the colour changed from colourless to pink. Then from the burette reading the amount of iron present in Mohr's salt solution was determined along

②

Procedure:-

1. 0.001 M NaBH_4 was prepared with distilled water in a 100 ml volumetric flask.
2. 0.0001 M AgNO_3 was prepared with distilled water in a 100 ml volumetric flask and taken 50 ml in an Erlenmeyer flask.
3. A magnetic stirrer bar was added in the flask and stirred it with a constant speed.
4. NaBH_4 solution was added drop wise in the flask with stirring.
5. After yellow coloration ~~a~~ excess NaBH_4 solution was added in the solution.

Observation:-

In 0.0001 M AgNO_3 solution, 0.001 M NaBH_4 ^{solⁿ} ~~was~~ will be added with dropper. After few drop of NaBH_4 around 7-8 drops it will turn yellow and after addition of more drops it will turn golden. Further, if more drops of NaBH_4 solⁿ will turn into blackish colour. If it will keep for around 12 hours long it will turn into agglomeration of silver nanoparticles.

① Experiment no:- 3

Date:- 02-11-23

Synthesis of Silver Nanoparticles Doped Organic Thin Film For Organic Transistors:-

AIM:- Synthesis of Silver Nanoparticles through reduction route.

Theory:- The definition of nanoparticles is an aggregate of atoms bonded together with a radius between 1 and 100 nm. There is a sudden shift of all properties of material when they just enter into the nanoscale. As material size reduces from centimeter (bulk) to nanometer scale, properties mostly decreases as much as six orders of magnitude to that at macro level. As nanostructures are having reduced dimensions, it leads to increase in surface energy via increase in surface area. If we consider bulk material the surface to volume ratio is low, whereas the surface to volume ratio is increasing enormously in case of nanomaterial. Conducting metals may show insulating behaviour in nanodimension. Hence one nanometer is a magical point on the dimensional scale.

Materials needed:-


Silver nitrate (AgNO_3); Sodium Borohydride (NaBH_4),
Distilled water.

⑥
Conclusion :-


The amount of Iron present in Mohr's salt solⁿ is 5.5428 g/L.

Date of Experiment :- 12-10-23

Date of Submission :- 2-11-23

Signature of Teacher :-  2/11/23

③



The quick complex and one step interaction of polydentate ligands with metals to yield stable complexes, is the principle of these complexometric titrations in which metal in solutions are titrated against that of polydentate ligands.

Apparatus →

- Pipette (25ml)
- Burette
- Beaker (250ml)
- Conical Flask

Reagents →

- Water sample
- EDTA solution
- Eriochrome Black-T indicator
- Buffer solution ($\text{NH}_4\text{Cl} + \text{NH}_4\text{OH}$)



Procedure :-

Determination of total hardness :

50 ml of water sample is taken in 250 ml conical flask. 1 ml of buffer solution is added to it. Now a pinch of EBT indicator is added. The solution becomes wine red. Now the solution

it is used to estimate Ca^{2+} & Mg^{2+} ions. Dissociation constant values indicate the EDTA behaves like a dicarboxylic acid. Two of its carboxyl group are strongly acidic the other two hydrogens are released during complex formation.

Tonisation of this complex depends on the pH of the solution. Hence in this titration pH of the sensitive indicators are used. In this estimation of Ca^{2+} and Mg^{2+} with EDTA, an azodye called Erichrome Black-T (EBT) is used as an indicator. This forms a metal indicator complex the stability is initially red. As the titration progress the metal ions form more stable complex with EDTA, hence the indicator anion goes to the solution and the colour changes from wine red to blue at the end point. Since the action of the indicator of formation of metal-EDTA complex is governed by pH hence it is kept constant by adding a suitable buffer ($\text{NH}_4\text{Cl} + \text{NH}_4\text{OH}$).



①

Experiment no → 4

Date → 09-11-23

Determination of Water Quality Parameter (Hardness of water)

AIM →

Estimation of calcium and magnesium hardness separately of a given tap water sample.

THEORY →

Hardness is defined as soap consuming capacity of water. Bicarbonates of Ca, Mg cause temporary hardness. However, chlorides, sulphates and carbonates of Ca & Mg cause permanent hardness. Hard water is not suitable for use as boiler feed water, because it leads to the formation of scales & slugs causing wastage of fuel, chocking of pipes, decreases efficiency & chances of explosion is also present. Hard water may also produce undesirable spots on the fabric due to its reaction with dyes when using in textile finishing. Thus softening of water is very important. In order to use any softening process type & extent of hardness must be known as prerequisite; complexometric titration is one of the best methods for hardness estimation.

Many metals form complexes with such reagents which contain appropriate ligands. EDTA is one such reagent, which forms complexes with metals. In the form of its disodium salt,

Conclusion

At first if 0.001 M NaBH_4 solution will added dropwise in 0.0001 M AgNO_3 solution; the colour would turn out as yellow colour.

Date of Experiment: 2-11-23

Date of submission:- 9-11-23

Signature of Teacher: J. Dutta

9/11/23

Experiment - 5

Date → 23-11-23

Determination of the concentration of the Electrolyte Through pH Measurement.

AIM:

Determination of the unknown strength of HCl solution by Standardised NaOH solution Using pH Metric Method.

THEORY:

Most of the chemical & biochemical processes are profoundly affected by the acidity and / or the alkalinity of the medium in which they take place. All acids dissociate in aqueous solution to yield H^+ ions. Some acids such as HCl, H_2SO_4 and HNO_3 are completely ionized in aqueous solution, whereas most of organic acids like $HCOOH$, CH_3COOH , ionized to a small extent only. The former are known as 'strong' acids, while the latter are classified as 'weak' acids. If we compare the ionization of 1(N) HCl and 1(N) CH_3COOH , then we will notice that, even their molar concentrations are same. HCl is ionized extensively and forms a strong acid solution whereas the CH_3COOH forms a weak acid solution containing lesser amount of H^+ ions. From these observations, we can conclude that the acidity of a solution does not depend upon the molarity of the acid, but upon the concentration of H^+ ions. We can follow from the above discussion that

⑥

$= 244 \text{ mg of } \text{CaCO}_3$

Therefore the total Hardness of water sample is $= 244 \text{ ppm}$.

Discussion

1. NH_4Cl and NH_4OH buffer solution is added to maintain the pH of the reaction medium at 10 which is a prerequisite for formation of metal - EBT complexes.
2. After addition of EDTA solution stable metal - EDTA complex is formed. Thus free EBT is responsible for the blue colouration of the solution after the titration is over.

Conclusion

Total Hardness of water sample $= 244 \text{ ppm}$.

Date of Experiment $\rightarrow 09-11-23$

Date of submission $\rightarrow 30-11-23$

Signature of teacher $\rightarrow Y. Dutta$

30.11.23

⑤

Result and Calculation

Table - I → Determination of Total Hardness :-

No of Observations	Volume of water sample taken (ml)	Burette reading		Volume of EDTA sol ⁿ consumed (ml)	Concordant reading (ml)	Hardness in ppm
		Initial	Final			
1	50	0.0	12.2	12.2	12.2	244
2		12.2	24.4	12.2		
3		24.4	36.8	12.4 <i>Sum 36.8 - 24.4 = 12.4</i>		

1000 ml of 1 M EDTA solution = 100 g of CaCO_3

1 ml of 0.01 M EDTA = $\frac{100 \times 0.01}{100}$ g of CaCO_3 .

12.2 ml 0.01 M EDTA = $\frac{100 \times 0.01 \times 12.2 \times 1000}{1000 \times 100}$

= 12.2 mg of CaCO_3 .

50 ml of Tap water sample contains 12.2 mg of CaCO_3 .

1000 ml " " " " " $\frac{12.2 \times 1000}{50}$

is titrated against 0.01(M) EDTA solution. At the end point the colour changes from wine red to blue. The titration is repeated three times and reading is noted.

Table → Titration of HCl using NaOH →


Serial No.	No. of Drops of NaOH solution	pH
1	0	2.96
2	3	3.02
3	6	3.09
4	9	3.13
5	12	3.26
6	15	3.34
7	18	6.12
8	21	9.23
9	24	9.84
10	27	10.11
11	30	10.36
12	33	10.52
13	36	10.70
14	39	10.87
15	42	11.07

Drop counting:

20 drops = 1 ml

1 " = $\frac{1}{20}$ ml

No of Drops of NaOH solⁿ Vs pH graph.



③
Whereas a_{H^+} is the activity of the H^+ ions. For dilute solutions, activity can be replaced by concentration, i.e.,


$$pH = -\log_{10}[H^+]$$

All pH meters have provision for standardizing the glass electrode in a buffer solution of known pH. This is necessary because different electrodes have different asymmetry potentials. One adjustment has been made, so that the registers correctly the known pH of the buffer solution, the instrument gives the pH of other solution without any calculation.

Measurement of pH is also employed to monitor the course of acid base titrations. The pH value of the solution at different stages of acid base neutralization is determined and is plotted against the volumes of the acid / alkali added. On adding a base to an acid, the pH rises slowly in the initial stages $pH = -\log_{10}[H^+]$; then it changes rapidly. However, the shape of the inflection point (i.e., where the pH changes abruptly) and symmetry of the curve on its two sides depends upon the ionizability of the acid and the base used and on the molarity of the acid and the acidity of the base.

if we know that about the concentration of H^+ ions. i.e. $[H^+]$, for a given solution the acidity of that solution can be expressed qualitatively.

Water is weakly ionized into H^+ ions and OH^- ions according to the following reaction:-



The ionic product of water, $K_w = [H^+][OH^-]$

For pure water $[H^+] = [OH^-] = 10^{-7}$

The H^+ ion concentration can be measured accurately using suitable methods. The acidic and alkaline nature of a solution can be expressed in the terms of concentration of H^+ ions. Due to some practical difficulties, the acidity or alkalinity of a solution is not expressed in terms of concentration of H^+ ions. Sorensen suggested the use of pH scale.

H^+	10^0	10^{-1}	10^{-2}	10^{-3}	10^{-4}	10^{-5}	10^{-6}	10^{-7}	10^{-8}	10^{-9}	10^{-10}	10^{-11}	10^{-12}	10^{-13}	10^{-14}
pH	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14

← Acidic → Neutral ← Alkaline →

pH scale is a logarithmic scale. It is defined as the -ve logarithmic of the H^+ ion concentration. The pH scale has values between 0-14. On this scale, pH of any solution is defined as:-

$$pH = -\log_{10} a_{H^+}$$

Aparatus :-

- Burette
- Pipette (10ml)
- Digital pH Meter & pH Electrode
- Plastic Beaker (100ml)

Reagents :-

- Unknown HCl solution
- Standardized NaOH ($\frac{N}{10}$) solution [supplied]

Procedure :-

1. Standardized NaOH ($\frac{N}{10}$) solution is supplied.
2. Unknown HCl solution was supplied.
3. The pH electrode was rinsed with de-ionized water.
4. 10 C.C HCl solution was pipetted out into the plastic beaker and water was added if necessary, so that both the electrodes were completely immersed within the solution. The electrode was joined with the Digital pH Meter and the pH was measured very carefully.
5. NaOH solution was added from a bottle drop-wise (approximately 2-3 drops)
6. The pH value of the solution was measured after addition of 2-3 drops of NaOH and the pH electrode was mildly stirred. The process was repeated until you had at least

6

Calculation V_1 = Volume of NaOH solution (result from graph V_2 = Volume of HCl = 10 ml $\times \frac{1}{20}$ ml) S_1 = Strength of NaOH = 0.11 (N) S_2 = Strength of HCl = ?

$$V_1 S_1 = V_2 S_2$$

$$S_2 = \frac{V_1 \times S_1}{V_2}$$

$$\Rightarrow S_2 = \frac{19.8 \times \frac{1}{20} \times 0.11}{10} = \frac{0.99 \times 0.11}{10} \text{ (N)}$$


$$\therefore S_2 = 0.01 \text{ (N)}$$

Conclusion:

The unknown strength of HCl is determined by standardized NaOH solution with the help of a digital pH meter. Hence the unknown strength of HCl is - 0.01 (N).

Date of Experiment : 23-11-23

Date of Submission : 30-12-23

Signature of Teacher: 

5
five points beyond the end point.

7. A curve was drawn by plotting the pH value against the no of drops of the titrant, the end point was found and the strength of the HCl solution was calculated.

Discussion :-

1. NaOH is a secondary standard solution, so it should be standardized with primary standard Oxalic acid solution using formula

$$V_1 S_1 = V_2 S_2$$

2. No indicator is used during the titration of HCl against NaOH solution.

3. This experiment is more accurate as it is performed by a digital pH meter. Thus we can minimize human error.

4. After plotting the graph (drops of NaOH solution along x-axis and pH along y-axis), we get the value of NaOH solution. So the strength of unknown HCl solution will be more accurate than other titration.

5. Titration is performed at room temperature (25° - 30° C)